

ACCESSORIES

INFUSION STANDS

- Stainless steel construction
- Rail mounted
- Height adjustable
- Various length turning arm

MOUNTED CABINETS

- Stainless steel and painted construction
- 6 models, 3 types of mounting
- Convenient for cleaning - without sharp edges
- Operative mounting to the rails or console profiles

EXAMINATION LIGHTS

- Power 25/ 35 W
- Voltage 12 V/220 V
- Light intensity in 1 m distance ~ 8000 lux
- Source setout range 1 m
- Railing mounted

RAILINGS FOR MOUNTED DEVICES

- 2 types of mounting: wall mounted or console profile mounted
- Dimensions for device mounting 25 x 10 mm
- Lifting capacity 40kg/m

BEDSIDE SCREENS

- Two types of holders:
 - telescopic
 - ceiling profile
- Wall rail or ceiling mounted
- Various length and height
- Easy cleaning, flame resistant material

Procedural tables

-
- Stainless steel and painted construction

- With drawer or without drawer

- Easily removable stainless steel trays
- Various dimensions

