

OTHER PRODUCTS

DENTAL LABORATORY EQUIPMENT

- dental technician's

workstation

-

Comfortable light source

-

Double dust filtration & collection

-

Dust collection system synchronous with micromotor

-

Low noise, < 59 dB (A)

-

Air blower with non-straining and auto-return hose

-

Gas armature

-

Stainless steel suction intake

-

Comfortable removable hand supports

-

Table measurements: 1200 x 630 x 850 mm

EDUCATIONAL EQUIPMENT

- dental simulation

equipment

-

-

-

Functional possibilities of Dental unit

-

All directions adjustable halogen operating lights

-

Universal control system for 3 instruments

-

Various phantom heads available

-

Universal suction system

-

All directions adjustable simulator

- Table measurements: 1200 x 630 x 850 mm

DENTAL CLINIC EQUIPMENT

- oil-less medical air

compressor

-

100% oil-less

-

0,01 μm filtration

-

Pressure regulator with monometer

-

One or two phases power 240/380 V 50Hz

-

Dew point up to -40°C

-

Low noise < 60db (A)

-

Stainless steel tank and dryer

-

Productivity 90 -230 l/min

